Agilent 6400 Series Triple Quadrupole LC/MS

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

G6475A - Control Panel

MANAGEMENT

Edit Instrument Delete Instrument Refresh All Lock Instrument **Configure Instrument** Create Acquisition Desktop Shortcuts Close Connection Copy to Clipboard Offline Worklist Editor Method Comparison Viewer Study Manager

Instruments and Locations Actions Selected Row LC/MS Acquisition Tools

1

MANAGEMENT

Create Edit Location Delete Location Refresh All

Instruments and Locations

Instruments

DESKTOP-9CKQSH3

Instruments Projects Administration

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
Role Membership

Password: Create Password

Confirm password: Confirm Password

Full name: Operator name in report

Email:

Contact Information:

☒ User must change password at next logon ← Unmark

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
Name (ID):

Description:


General Group Membership **Role Membership**

Member of:

	Name	Description
☒	Everything	All privileges
☐	System Administrator	Manage users and security settings
☐	Instrument Administrator	Manage instruments and locations
☐	Project Administrator	Manage projects and project groups
☐	Instrument User	View and run instruments
☐	Technician	Laboratory technician
☐	Chemist	Analytical chemist

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Create User Edit User Remove User

Administration

- My Settings
- Local Configuration
- System Configuration
- Security Policy
- Users**
- Groups
- Roles
- System Activity Log
- Links
- Licenses
- Instrument Controllers
- Diagnostics
- Instrument Options
- Administrative Reports

Users

Name	Description	Full name
User		User
admin		admin

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

MS Checktune Report - G6475A


Instrument Information

Model	G6475A	Checktune Date	2022-11-08T15:08:34-08:00
Serial Number	SG2222S002	SW/FW Version	3.0.1463/8.1.38
Ion Source	AJS ESI	Ionization Mode	ESI
Last Autotune Date	2022-11-08T11:26:26-08:00	Last Tuned By	SYSTEM (SYSTEM)
Vacuum Pressure	1.87E+0 [R] (Torr); 3.75E-5 [H] (Torr)	Overall Result	Passed

Positive Ion Mode

MS1 Peak Width Unit, Scan Speed Normal	Result	Passed
MS2 Peak Width Unit, Scan Speed Normal	Result	Passed
MS1 Peak Width Narrow, Scan Speed Normal	Result	Passed
MS2 Peak Width Narrow, Scan Speed Normal	Result	Passed
MS1 Peak Width Wide, Scan Speed Normal	Result	Passed
MS2 Peak Width Wide, Scan Speed Normal	Result	Passed
MS1 Peak Width Widest, Scan Speed Normal	Result	Passed
MS2 Peak Width Widest, Scan Speed Normal	Result	Passed
MS2 Scan Speed Fast	Result	Passed
MS2 Scan Speed Ultra	Result	Passed
MS1 Lag Factor	Result	Passed
MS2 Lag Factor	Result	Passed
Gain	Result	Passed

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
The screenshot shows the 'MS Checktune Report - Ultivo' window. The report includes the following sections:

- Instrument Information:**
 - Model: 004558 Ultivo
 - Serial Number: 50118452203
 - Ion Source: AJS ESI
 - Last Autotune Date: 2023-01-17T15:47:31-08:00
 - Overall Result: Passed
 - Checktune Date: 2023-01-17T16:25:34-08:00
 - SW/PS Version: 3.0.14675.1.38
 - Ionization Mode: ESI
 - Last Tuned By: SYSTEM (SYSTEM)
- Vacuum And Temperature:**
 - Rough Vac (Torr): 1.88E+0
 - Turbol1 Speed (%): 100.0
 - MS2 Heater (°C): 65
 - High Vac (Torr): 2.62E-5
 - MS1 Heater (°C): 65
- Positive Results:**
 - Components:**
 - Ion Source Settings:**
 - Gas Temperature (°C): 300
 - Nebulizer (psi): 15
 - Sheath Gas Flow (L/min): 11.0
 - Nicote Voltage (V): 1500
 - Optics Settings:**
 - Fragmentor (V): 120
 - Cyclone Outer DC (V): 3
 - Cyclone Exit DC (V): 3
 - Cyclone Inner RF (V): 400
 - Gas Flow (L/min): 7.0
 - Sheath Gas Temperature (°C): 200
 - Capillary Voltage (V): 4000
 - Skimmer (V): 10
 - Cyclone Accelerator (V): 0
 - Cyclone Exit Lens (V): 2
 - Cyclone Outer RF (V): 300

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

The screenshot shows the 'Method Editor' window. The 'TQ' tab is selected in the top toolbar. The 'Autotune' button is highlighted in the 'Tune' section. The 'Acquisition Parameters' section shows the following settings:

- Compound group:** sulfachloropyridazine, sulfadimethoxine, sulfamethazine, sulfamethizole
- Estimated cycle time (ms/cycle):** 28

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MS Autotune Report - G6475A

Instrument Information

Model	G6475A	Autotune Date	2022-11-08T11:26:26-08:00
Serial Number	SG22225002	SW/FW Version	3.0.1463/8.1.38
Ion Source	AJS ESI	Ionization Mode	ESI
Last Tuned By	SYSTEM (SYSTEM)	Overall Result	Passed

Vacuum And Temperature

Rough Vac (Torr)	1.91E+0	High Vac (Torr)	3.81E-5
Turbo1 Speed (%)	100.0	Turbo2 Speed (%)	100.0
MS1 Heater (°C)	100	MS2 Heater (°C)	100

Positive Ion Mode


Detector Settings

HED Standard (kV)	-10	EMV Standard (V)	1201
HED Large Molecule (kV)	-20	EMV Large Molecule (V)	1117

MS1 Peak Width Unit, Scan Speed Normal

Calibration Information		Result	
Mass Calibration Results		Passed	
Theoretical (m/z)	Mass Delta(m/z)	Width Delta(m/z)	Abundance
58.07	0.00	0.00	9.35E+5
118.09	0.01	-0.01	2.01E+6
322.05	0.00	0.00	4.87E+5
622.03	0.01	-0.01	5.68E+5
922.01	0.00	0.00	4.79E+5
1221.99	0.00	0.00	2.82E+5
1521.97	0.01	0.00	1.42E+5
1821.95	0.00	0.00	5.92E+4
2121.93	0.00	-0.01	4.32E+4

MS2 Peak Width Unit, Scan Speed Normal

Checktune vs Autotune



Checktune Weekly

- Routine use only requires resolution and mass axis verification for polarity and resolution modes.
- Checktune checks all three resolution modes.
- Checktune checks current tune file, but does not change it. No need to save Tune file.


Autotune Monthly/Quarterly

- Autotune in both polarities to optimize ion transmission and update EM voltage.

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Project: G6475A Training_Nov3, Instrument: G6475A - Acquisition

Home | Training | Method | Single Sample | Status | Save | Delete | Lock Layout | Reset | Instrument Status | Chromatogram Plot | Worklist | Method Editor | Sample Run | Spectrum | Method Optimizer | dMRM Method Split | Analytical Column Setup | Device Config | System Activity Log | Tools

For Help, Press F1 | No worklist loaded. | D:\Projects\...Basic.m | D:\ Free space: 277.7 GB

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

The screenshot displays the 'Instrument Status' window with five modules: Multisampler, Binary Pump, DAD, Column Comp., and TQ. Each module has a status indicator (Idle or Tune) and an EMF checkmark. A red line highlights the 'Instrument Tune' button at the bottom right, which is currently set to 'On'. Below the main window, three context menus are shown, each corresponding to a module's quick access options:

- Multisampler Context Menu:**
 - Control...
 - Error Method
 - Identify Device
 - Home All
 - Reset Injector
 - Switch to Mainpass Left
 - Switch off Tray Illumination
 - Auto-clean...
 - Prime Solvents...
 - Start Purge
 - Modify
 - Assign Wellplates
- Binary Pump Context Menu:**
 - Control...
 - Error Method
 - Identify Device
 - Switch On
 - Switch Solvent Selection Valve A
 - Switch Solvent Selection Valve B
 - Bottle Fillings...
 - Prepare Pump...
 - Seal Wash Prime
- Column Comp. Context Menu:**
 - Control...
 - Error Method
 - Identify Device
 - Switch Off
 - Column Assignment...

At the bottom, a detailed 'Instrument Tune' menu is shown, with the 'On' option selected. The 'Waste' option is also highlighted in the 'MS' section.

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Properties DA Multisampler M

Method

- Acquisition
- Source
- Chromatograms
- Timetable
- MRM database browser
- Convert to dMRM**
- Update dMRM method
- Advanced override settings

Acquisition Parameters

cursor m/z	MS1 res	Product m/z	MS2 res	RT (min)	RT window (min)	Fragmentor (V)	CAV (V)	CE (V)	Average dwell (ms)	Polarity
.1	Unit	152.2	Unit	0	1	105	2	12	17.36	Positive
.1	Unit	137.2	Unit	0	1	105	2	24	17.36	Positive
.1	Unit	174.1	Unit	0	1	125	3	16	17.36	Positive
.1	Unit	68.0	Unit	0	1	125	3	40	17.36	Positive
.1	Unit	165.1	Unit	0	1	80	2	20	17.36	Positive
.1	Unit	123.1	Unit	0	1	80	2	30	17.36	Positive
.1	Unit	169.1	Unit	0	1	105	2	32	17.36	Positive
.1	Unit	97.0	Unit	0	1	105	2	40	17.36	Positive
.0	Unit	198.8	Unit	0	1	70	5	0	17.36	Positive
.0	Unit	125.0	Unit	0	1	70	5	16	17.36	Positive
.1	Unit	69.0	Unit	0	1	115	2	16	17.36	Positive

Method Editor

Properties DA Multisampler Multisampler Pretreatment Binary Pump DAD Column Comp. TQ

Method

- Acquisition
- Source
- Chromatograms
- Timetable
- MRM database browser
- Update dMRM method**
- Advanced override settings
- Optimization Settings
 - Parameters
- Tune
 - Autotune

Select MassHunter TQ data file or Quant report folder

D:\Projects\G6475A Training_Nov3\Data\Pesticides_TestMix_MRM.d

Add additional compounds/transitions

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
The image shows two screenshots of the Agilent Worklist Run Parameters dialog box. The first screenshot shows the 'Run Parameters' tab with the 'Post-worklist' checkbox highlighted. The second screenshot shows the 'Data File Settings' tab with the 'Root file Folder' highlighted.

Worklist Run Parameters

Run Parameters | Data File Settings | iReflex | Additional Parameters

Operator Information

Operator name: SYSTEM (SYSTEM)

Run Information

Part of Method to run: Both Acquisition and DA (dropdown) ☒ Synchronous DA

Run Type: Standard Start (dropdown) ☒ Stop worklist on DA error

Method Paths

Method: D:\Projects\G6475A Training_Nov3\Methods (dropdown)

Override DA: D:\Projects\G6475A Training_Nov3\Methods (dropdown)

Scripts

☐ Pre-worklist (dropdown)

☐ Post-worklist (dropdown)

☐ Acquisition clean-up (dropdown)

☐ Post-analysis (dropdown)

Disk space Information (GB)

277 Available Disk space 200 Free Disk space threshold

Run Settings

☒ Overlapped injection Wait Time for Ready 10 (min)

Worklist Run Parameters

Run Parameters | **Data File Settings** | iReflex | Additional Parameters

Root Folder

Root file Folder: D:\Projects\TQ\Data (dropdown)

Sub-folders

Sub-folder 1: <Empty> (dropdown) Sub-folder 2: <Empty> (dropdown) Sub-folder 3: <Empty> (dropdown)

File Naming

Part 1: User Text (dropdown) Part 2: Counter[0001] (dropdown) Part 3: <Empty> (dropdown)

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





Table 1 LC parameters for sulfa drug mix

Parameter	Value for all instruments
PUMP	
• Flowrate	600 µL/min
• Solvent A	Water with 0.1% Formic Acid
• Solvent B	Acetonitrile with 0.1% Formic Acid
• Gradient (minutes - %B)	Initial Conditions: 90% Channel A and 10% Channel B 0 minutes - 10% B 5.0 minutes - 90% B
• Stop Time	5 minutes
• Post Time	3 minutes
INJECTOR	
• Inj. Vol.	1 µL
• Injection	Standard
• Draw Position	3.0 mm
UV DETECTOR	
• Ch A	272 nm (100 nm BW on DAD)
• REF A (DAD only)	360 nm (100 nm BW)
COL THERM	
• Temp	40° C

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Table 3 MS parameters for sulfa drug mix

Parameter	Value	
• Inlet	ESI (positive polarity)	ESI (positive polarity) with Agilent Jet Stream Technology
• Gas Temp	350 °C	350 °C
• Scan Type	MS2Scan	MS2Scan
• Nebulizer	50 psi	35 psi (nitrogen)
• Dry Gas	12 L/min	10 L/min
• Range	100 to 400	100 to 400
• Sheath Gas Temperature	not applicable	400 °C
• Sheath Gas Flow	not applicable	12 L/min
• Nozzle Voltage	not applicable	0 V
• Capillary voltage positive	4000 V	4000 V
• Delta EMV pos	400 V	200 V

កម្រិត Scan type = Product Ion Scan

Acquisition Parameters

Ion source: AJS ESI

Stop time: ☐ As pump/No limit ☐ Limit (min): 1

☒ Time filter window (min): 0.07

Time Segments

Start time (min)	Scan type	Detector Gain Factor (+)	Detector Gain Factor (-)	Stored?
0	Product Ion	1	1	☒

Estimated cycle time (ms/cycle): 1047

Estimated max scan speed (Da/s): 1041

Data storage: Centroid

Threshold: 0

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Acquisition Parameters

Ion source: AJS ESI

Stop time: ☐ As pump/No limit ☐ Limit (min): 1

☒ Time filter window (min): 0.07


Time Segments

Start time (min)	Scan type	Detector Gain Factor (+)	Detector Gain Factor (-)	Stored?
0	MRM	1	1	☒

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


☐ Determine parameter range from method source setpoints

☒ Gas Temperature (°C)

Pre-wait (min) 30 Replicate(s) 1 Step Wait (min) 20

From 300 To 350 Step Size 10

☒ Sheath Gas Temperature (°C)

Pre-wait (min) 30 Replicate(s) 1 Step Wait (min) 20

From 350 To 400 Step Size 10

☒ Capillary Voltage (V)

Pre-wait (min) 0 Replicate(s) 1 Step Wait (min) 0

From 3600 To 4400 Step Size 200

☒ Nebulizer (psi)

Pre-wait (min) 0 Replicate(s) 1 Step Wait (min) 0

From 33 To 37 Step Size 1

☒ Gas Flow (L/min)

Pre-wait (min) 0 Replicate(s) 1 Step Wait (min) 0

From 9 To 11 Step Size 0.5

☒ Sheath Gas Flow (L/min)

Pre-wait (min) 0 Replicate(s) 1 Step Wait (min) 0

From 9.5 To 12 Step Size 0.5

☐ Nozzle Voltage (V)


Pre-wait (min) 0 Replicate(s) 1 Step Wait (min) 0

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
Type	Pre-Wait (min)	Replicates	Step Wait (min)	From	To	Step Size
Capillary	0	1	0	1000	6000	500
Nozzle Voltage	0	1	0	0	2000	500
Drying Gas Temp	20	2	10	150	350	50
Drying Gas Flow	0	1	0	Non-funnel 10	Non-funnel 13	1
				Funnel 17	Funnel 20	
Sheath Gas Temp	20	2	10	150	400	50
Sheath Gas Flow	0	1	0	10	12	1
Nebulizer	0	1	0	25	60	5

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



PI_01.d

Survey_01.d

Agilent MassHunter Qualitative Analysis software

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

Sample	Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time
SulfamixMR...	SulfamixMR...					

Quantifier	Name	Δ	TS	Transition	Scan	Type	Precursor Ion	Product Ion	RT	Ion Polarity
Compound_1	1	311.0	->	155...	MRM	Target	311.0	155.7	2.213	Positive
Compound_2	1	285.0	->	155...	MRM	Target	271.0	155.7	0.833	Positive
Compound_3	1	279.0	->	185...	MRM	Target	279.0	185.7	1.190	Positive
Compound_4	1	271.0	->	155...	MRM	Target	285.0	155.8	0.445	Positive

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



Method table

Time Segment: < All > Compound: AA Reset Table View Level Name Prefix: # of Levels: 3 Create Levels

Sample						
Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time	
SulfamixMR...	SulfamixMR...					

Quantifier							
Name	TS	Transition	Scan	Type	Units	Dil. High Conc.	Dil. Pattern
AA	1	311.0 -> 155...	MRM	Target	ng/ml	100.0000	1:2
BB	1	285.0 -> 155...	MRM	Target	ng/ml	100.0000	1:2
CC	1	279.0 -> 185...	MRM	Target	ng/ml	100.0000	1:2
DD	1	271.0 -> 155...	MRM	ISTD	ng/ml	100.0000	1:2

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Method table

Time Segment: < All > Compound: AA Reset Table View Level Name Prefix: # of Levels: 3 Create Levels

Sample						
Name	Data File	Type	Level	Acq. Method File	Acq. Date-Time	

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Overview

Startup and Shutdown

Calibrant Delivery
System and Divert
Valve

Sources

Nebulizer

Corona Needle

Ion Transfer Capillary

High Pressure Ion
Funnel

Desolvation Assembly
and Ion Optics

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



Error

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The screenshot shows the Agilent MassHunter Workstation Data Acquisition interface. In the top right corner, the 'System Activity Log' icon is highlighted with a red box. A red arrow points from this icon to the 'Agilent MassHunter Activity Log Viewer' window. The window displays a table of activity logs with columns for Date/Time, User, and Description.

Date/Time	User	Description
2021-06-28 16:47:07:00	SYSTEM (SYSTEM)	Activity Log Viewer: User 'SYSTEM' has logged in.
2021-06-28 16:19:22:07:00	SYSTEM (SYSTEM)	Acq Console: Context is switched to 'Tune'.
2021-06-28 16:19:22:07:00	SYSTEM (SYSTEM)	Acq Console: User has loaded layout 'Status(sys)'.
2021-06-28 16:19:11:07:00	SYSTEM (SYSTEM)	Acq Console: User 'SYSTEM' has logged in.
2021-06-28 16:19:09:07:00	SYSTEM (SYSTEM)	Acquisition Engine: Method 'C:\Projects\Test01_0412\Methods\default.m (Version: NA)' loaded on instrument.
2021-06-28 16:19:07:07:00	SYSTEM (SYSTEM)	Acq Console: User has loaded layout 'Method(sys)'.
2021-06-28 16:17:58:07:00	SYSTEM (SYSTEM)	Offline Method Editor: User 'SYSTEM' has logged out.
2021-06-28 15:50:59:07:00	SYSTEM (SYSTEM)	Offline Method Editor: User 'SYSTEM' has logged in.

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